



Albuterol Sulfate
Inhalation Solution 0.083%*

Artwork Type: **PIL**

Artwork Code: **PHPI0033**

Void Code:

Dimension: **93x215 mm**

Void A/W Reason:

Country: **ANDA-US**

Language: **ENGLISH**

Mfg. Location: **SPML**

Specification/Type of Paper:

SUPER FINE 41 GSM ITC PAPER

Folding:

93---0---93 mm

215---3---26.8 mm

Special Req.:

Remark (if any):

Prepared by: **SAPNA**

Checked by:

Approved by:

Approved by RA:

APPROVAL HISTORY ATTACHED

No. of Colors: 1

BLACK



PHPI0033

Albuterol Sulfate
Inhalation Solution

Patient's Instructions for Use


Albuterol (al-BUE-ter-ol) Sulfate Inhalation Solution 0.083%*

(*Potency expressed as albuterol)


Note: This is a unit-dose vial. No dilution is required.

Read complete instructions carefully before using.

1. Remove vial from the foil pouch.
2. Twist the cap completely off the vial and squeeze the contents into the nebulizer reservoir (Figure 1).



3. Connect the nebulizer reservoir to the mouthpiece or face mask (Figure 2).



Front side 93 mm


Size: 93x215 mm

Folding

93---0---93 mm

215---3---26.8 mm

4. Connect the nebulizer to the compressor.
5. Sit in a comfortable, upright position; place the mouthpiece in your mouth (Figure 3) (or put on the face mask); and turn on the compressor.



6. **Breathe as calmly, deeply and evenly** as possible until no more mist is formed in the nebulizer chamber (about 5 to 15 minutes). At this point, the treatment is finished.
7. Clean the nebulizer (see manufacturer's instructions).

Note: Use only as directed by your physician. More frequent administration or higher doses are not recommended.

Store Albuterol Sulfate Inhalation Solution 0.083% at 20° to 25° C (68° to 77°F); excursions permitted between 15° to 30°C (59° to 86°F). Store in pouch until time of use.

ADDITIONAL INSTRUCTIONS: _____

Distributed by:
Sun Pharmaceutical Industries, Inc.
Cranbury, NJ 08512



Manufactured by:
Sun Pharmaceutical Medicare Limited
Baska Ujeti Road, Ujeti
Halol-389350, Gujarat, India.

PHP10033
ISS. 01/2021

Back side 93 mm

Folding
93---0---93 mm
215---3---26.8 mm

Size: 93x215 mm